

# AWT AutoClave Water is a high-end water purification system to use with dental autoclaves

**AWT AutoClave Water delivers water to even the most demanding autoclaves such as DAC and other systems. It has a high reliability and a long running time between changing of consumables.**

All components is of high quality to meet the standards of Adept Water Technologies products.




Fig. 1: AWT AutoClave Water


### AWT AutoClave Water systems:

	AWT AutoClave Water Classic	AWT AutoClave Water Smart
Dimensions: L x W x H in mm	420 x 255 x 435	420 x 255 x 435
Weight kg	13	15
Connection water Inlet / Outlet, JG mm	8 / 8	8 / 8
Drain JG mm	8	8
Temperature in / max C°	2 / 40	2 / 40
Input pressure max / min bar	6 / 2,5	2,5 / 1
Max water storage l	7	7
Max production LPD	150	150
Max TDS input ppm	2000	2000
Power connection V	None	110 / 220
Water Quality out µS/cm	<1	<1
Order no.	140028	140029

**Fig. 2:** AWT AutoClave Water Classic – is only run by the inlet water pressure. Operate with an inlet water pressure 2.5 bar or higher.



**Fig. 2:** AWT AutoClave Water Smart – has a booster bump installed. Operate with an inlet water pressure 1 bar or higher.



### AWT AutoClave Water accessories:

AWT AutoClave Water conductivity meter Order no. 140030	InLine conductivity /TDS meter for measuring the conductivity / TDS after the AWT AutoClave Water system
--	--

Always consult your local dealer before connecting the system to an autoclave.

Adept Water Technologies are not responsible for any failures or breakdown of autoclaves because of the water coming from the AWT AutoClave Water system.