

Water prefilters

Sometimes the tap water is of such a poor quality it affects the performance of your BacTerminator® Dental. In those cases we recommend to prefilter

the water before using it. Adept Water Technologies therefore has a range of smart prefilters to be used in case of low quality tap water.



Fig. 1: Point of entry prefilter HF25/HF35

Fig. 2: Sediment filter 5 micron

Fig. 3: Carbon filter

Fig. 4: UF filter

Data for prefilters:

Type:	Point of entry prefilter HF35	Point of entry prefilter HF25	Sediment filter 5 micron	Carbon filter	UF filter
Capacity l	47.000	37.000	2.000	2.000	2.000
Flow rate lpm	6,3	5,7	2	2	2
Max temperature °C	37,8	37,8	38	38	38
Max Pressure bar	8,5	8,5	4	4	4
Inlet/Outlet	3/8" NPT	3/8" NPT	3/8" Tube	3/8" Tube	3/8" Tube
Size H x W x D, cm	54 x 11 x 11	30 x 8 x 8	25 x 5 x 5	25 x 5 x 5	25 x 5 x 5
Description	Pre-filter for larger volumes of incoming water	Pre-filter for larger volumes of incoming water	Pre-filter before the BacTerminator® Dental	Pre-filter before the BacTerminator® Dental	Pre-filter before the BacTerminator® Dental
Order no.	140033	140032	140034	140035	140036

Capacity and flow rate are only indicators. These products are only filters, not water purifiers. Filters are NSF Standard, and are approved as Medical

Device accessories for the BacTerminator® Dental. Local guidelines and legislation regarding drinking water, may apply.