

## Bracket systems for BacTerminator® Dental series

If it is inconvenient to have a disinfection unit on the floor next to the dental unit, brackets are available as accessories to the BacTerminator® Dental product series. The Bracket system makes it possible to mount the BacTerminator® Dental unit on a wall.

### Wall Bracket for BacTerminator® Dental Compact

Part no.: 160107

H x W x D: 160 mm x 200 mm x 170 mm

Weight: 1.3 kg



### Wall Bracket for BacTerminator® Dental

Part no.: 160083

H x W x D: 450 mm x 370 mm x 180 mm

Weight: 7 kg



### Double Wall Bracket

Part no.: 160084

H x W x D: 450 mm x 370 mm x 350 mm

Weight: 11,5 kg

